

CHARACTERISTICS OF DEIONISED WATER

(A) DEFINITIONS OF PURIFIED WATER


Table 5 provides a definition of purified water quality in terms of the ionic impurities only; it lists three grades of deionised water - ultra-pure, pure and purified water - each with a specified range of resistivity/conductivity. Table 6 is a conversion chart for conductivity and resistivity; one is the reciprocal of the other.

(B) RESISTIVITY AND TEMPERATURE

The resistivity and conductivity of water are temperature-dependent. If the temperature rises by 1°C, the conductivity of mains water will increase by approximately 2%, whereas that of ultra-pure water will increase by up to 6%. It is standard practise to correct conductivity and resistivity to 25°C, a process which is carried out automatically by modern conductivity meters. Figure 22 plots the resistivity of ultra-pure water against temperature.

(C) pH MEASUREMENTS IN DEIONISED WATER

The pH of water is a measure of its acidity or alkalinity, and is defined as: $\text{pH} = -\log_{10} [\text{H}^+]$, where $[\text{H}^+]$ is the molar concentration of hydrogen ions. pH values are defined on a scale from 1 to 14, with ultra-pure water having a neutral pH of 7.0.



The pH of ultra-pure water is difficult to measure. Not only does high-purity water rapidly pick up contaminants - such as carbon dioxide (CO_2) - that affect its pH, but it also has a low conductance that can affect the accuracy of pH meters. For instance, absorption of just a few ppm of CO_2 can cause the pH of ultra-pure water to drop to 4.5, although the water is still of essentially high quality.

The most accurate estimation of the pH of ultra-pure water is obtained by measuring its resistivity; for a given resistivity, the pH must lie between certain limits. For example, if the resistivity is 10.0 $\text{M}\Omega\cdot\text{cm}$, the pH must lie between 6.6 and 7.6. The relationship between the resistivity and pH of high-purity water is shown in Figure 23.

Table 5: Definition of purified water quality (in terms of ionic impurities only)

MΩ .cm @ 25°C	Ω .cm @ 25°C	μS.cm @ 25°C	Total dissolved solids in parts per million	
Electrical Resistivity		Electrical Conductivity	TDS in ppm (approx.)	
18	18,000,000	0.0555	-	ULTRA-PURE WATER
10	10,000,000	0.1	-	
5	5,000,000	0.2	-	PURE WATER
2	2,000,000	0.5	-	
1	1,000,000	1.0	0.7	
0.5	500,000	2	1.4	PURIFIED WATER
0.2	200,000	5	3.5	
0.1	100,000	10	0.7	
0.05	50,000	20	14	
0.02	20,000	50	35	

Quality:	ULTRA-PURE WATER	PURE WATER	PURIFIED WATER
Typical Resistivity: or	10 - 18 MΩ.cm	1 - 10 MΩ.cm	1 - 0.02 MΩ.cm
Conductivity:	0.1 - 0.0555 μS/cm	1.0 - 0.1 μS/cm	1 - 50 μS/cm
Produced by:	Polishing mixed-bed system e.g. nuclear grade resins.	Strongly basic mixed-bed system.	Weakly basic mixed - bed system.

Table 6: Conductivity versus Resistivity Conversion Chart




Figure 22: Relationship between the resistivity of ultra-pure water and temperature




Figure 23: Electrical resistivity versus pH of deionised water

